Automatic Veterinary Syringes

2.0 cc ThaMa 253

Fix Dosage

(with 3 different pistons) 0.5cc, 1.0cc, 2.0cc With vial holder

Operating Instructions:

This Automatic Syringe ThaMa 223 is a precision, reliable and accurate All-Metal-Fixed-Dosage instrument, designed for mass treatment of poultry. All the components are made of the finest material, oil proof, non-corrosive; the valve is made of stainless steel ball. The piston slides smoothly along the metal cylinder.

It is supplied with 3 pistons of various sizes: 0.5cc, 1.0cc, 2.0cc.

All parts can be sterilized up to 125°C.

Before Use:

- 1. It's recommended to sterilize the syringe before each use.
- 2. Check that all screwed parts are tightened.
- 3. Make sure that the valve springs and washers are properly positioned into the base.

To set the dose:

- 1. In the box find the key-pin (20).
- 2. Open the body-cover (13) by gripping it with your fingers and turn.
- 3. Press the piston forward to its utmost position, and insert the key-pin into the hole on the piston.
- 4. Grip the piston and unscrew it, and change to the one you need.
- 5. Tight the new piston gently with the key-pin.
- 6. Take out the key-pin from the piston.
- 7. Put a drop of caster-oil on the piston's O-ring.
- 8. Close the body-cover and tight it.

Using the Vial Holder:

- 1. Remove the black plastic clip (42).
- 2. Insert the vial (43) to the vaccine bottle and push it down.
- 3. Return the plastic clip back connect the vial holder to the syringe.

Pay attention that the black plastic clip is properly set!

After use:

Care and Maintenance:

- 1. Dismantle the vial from the vial holder.
- 2. Pump the syringe in clear water, 6-10 times, to remove any residue of the material from the body, needles, tube and valves.
- 3. Open the body cover (13) and wash all parts in clear water.
- 4. Open the luer lock needle adapter (02) and the tubing valve housing (03) wash them in clear water.
- 5. Oil the glass cylinder after cleaning.
- 6. Assemble the parts together.

Note: Pay attention, and take good care not to lose the small parts, like valves, valves springs, and valves plastic washers.

Note: It is highly recommended to oil the glass cylinder after cleaning and before assemble the parts together

2.0 cc ThaMa 253 Fix Dosage

(with 3 different pistons) 0.5cc, 1.0cc, 2.0cc

With vial holder

Recommended:

Use the original box for each syringe to keep it useful for long time. The automatic syringe is made of the best materials, is oil-proof, non corrosive. All parts can be sterilized up to 125°C. The glass barrel is made of clear and tempered glass, exactly graduated, enabling correct setting and continuous control during operation.

253-01	Syringe Body
253-02	Luer lock needle adapter
253-05	Valve Spring
253-06	Valve stainless ball
253-07	Valve Plastic Washer
253-10	Piston O-Ring
253-11	Piston Shaft
253-13	Body Cover
253-14	Handle Spring
253-16	Finger Grip Handle
253-17	Upper Handle
253-18	Spring Housing
253-42	Plastic Clip
253-43	Vial Holder – Body
253-99	Spare Parts-Kit (include: 253-
05*2, 253-06	i. 253-07*2. 253-10*5)

