Automatic Veterinary Syringes

0.5 cc ThaMa 223 Fix Dosage

(with 5 different pistons) 0.1cc, 0.2cc, 0.25cc, 0.3cc, 0.5cc Using vial holder as feeding bottle



Operating Instructions:

This Automatic Syringe ThaMa 223 is a precision, reliable and accurate All-Metal-Fixed-Dosage instrument, designed for mass treatment of poultry. All the components are made of the finest material, oil proof, non-corrosive; the valve is made of stainless steel ball. The piston slides smoothly along the metal cylinder.

It is supplied with 5 pistons of various sizes: 0.1cc, 0.2cc, 0.25cc, 0.3cc, 0.5cc.

All parts can be sterilized up to 125°C.

Before Use:

- 1. It's recommended to sterilize the syringe before each use.
- 2. Check that all screwed parts are tightened.
- 3. Make sure that the valve springs and washers are properly positioned into the base.

To set the dose:

- 1. In the box find the key-pin (20).
- 2. Open the body-cover (13) by gripping it with your fingers and turn.
- Press the piston forward to its utmost position, and insert the key-pin into the hole on the piston.
- 4. Grip the piston and unscrew it, and change to the one you need.
- 5. Tight the new piston gently with the key-pin.
- 6. Take out the key-pin from the piston.

7. Put a drop of caster-oil on the piston's O-ring.

8. Close the body-cover and tight it.

Choose the injection system you need: Vial Holder:

- 1. Remove the black plastic clip (42).
- 2. Insert the vial (43) to the vaccine bottle and push it down.
- 3. Return the plastic clip back connect the vial holder to the syringe.

Pay attention that the black plastic clip is properly set!

Feeding Bottle:

- 1. Connect the feeding bottle to the syringe.
- 2. To fill the feeding bottle, just remove the "rear red stopper".
- 3. Fill the feeding bottle with the vaccine needed.
- 4. Now start working.
- 5. To refill the feeding bottle, just remove the rear red stopper.
- 6. For cleaning: fill with water and pump.

After use:

Care and Maintenance:

- Dismantle the vial from the vial holder.
- 2. Pump the syringe in clear water, 6-10 times, to remove any residue of the material from the body, needles, tube and valves.
- 3. Open the body cover (13) and wash all parts in clear water.
- Open the luer lock needle adapter (02) and the tubing valve housing (03) - wash them in clear water.
- 5. Oil the glass cylinder after cleaning.
- 6. Assemble the parts together.

Note: Pay attention, and take good care not to loose the small parts, like valves, valves springs, and valves plastic washers.

Note: It is highly recommended to oil the glass cylinder after cleaning and before assemble the parts together.

0.5 cc ThaMa 223 Fix Dosage

(with 5 different pistons) 0.1cc, 0.2cc, 0.25cc, 0.3cc, 0.5cc Using vial holder as feeding bottle

Recommended:

Use the original box for each syringe to keep it useful for long time. The automatic syringe is made of the best materials, is oil-proof, non corrosive. All parts can be sterilized up to 125°C. The glass barrel is made of clear and tempered glass, exactly graduated, enabling correct setting and continuous control during operation.

223-01	Syringe	Body
--------	---------	------

223-02 Luer lock needle adapter

223-05 Valve Spring

223-06 Valve stainless ball

223-07 Valve Plastic Washer

223-10 Piston O-Ring

223-11 Piston Shaft

223-13 Body Cover

223-15 Handle Spring

223-16 Finger Grip Handle

223-17 Upper Handle

223-18 Spring Housing

223-29 NBR valve

223-71 Valve fiber washer

Feeding bottle parts

223-42 Plastic Clip

223-43 Vial Holder – Body

223-44 Rear Feeding Stopper

223-45 Feeding Bottle – Body

223-46 Front Feeding Plug

223-47 Feeding Valve Housing

223-99 Spare Parts-Kit (include: 223-05*2, 223-06*2, 223-07*1, 223-10*4, 223-29*1, 223-71*2)

223-20 Key-pin

223-21 Piston 0.1cc

223-22 Piston 0.2cc

223-23 Piston 0.3cc

223-25 Piston 0.25cc

223-24 Piston 0.5cc

